

**NEW**

## READY DNA Extraction

### BEXS 12

Cost-Effective, Automatic DNA Extraction With  
Our New Ready DNA Bead Kit and BEXS 12

- 1 – 12 samples per run
- Standard protocol on a CF protocol card
- For 200µl of whole blood
- Run time approx. 50 minutes
- Yield: approx. 4 – 8µg of DNA
- Reagents can be connected directly
- Cost-efficient as without Cartridge System
- CE certified and optimized for use in molecular HLA and blood group detection



**inno-TRAIN**  
Diagnostik GmbH

QMS 02.08



## READY DNA Extraction

### DNA Extraction With Our New Ready DNA Bead Kit Is Also Suitable For Manual Handling

- Same bead technique as for the BEXS 12
- No centrifugation steps necessary
- Special Ready DNA Bead Stand for magnetic separation
- For 250µl of whole blood
- Extraction time approx. 30 minutes
- Yield: approx. 5 – 10µg of DNA
- CE certified and optimized for use in molecular HLA and blood group detection



#### Article:

GX01050	BEXS 12 Bead Extraction System	
GX01052	BEXS 12 Ready DNA Bead Kit	for 180 Extractions
GX01053	Ready DNA Bead Kit manual	for 50 Extractions
GX01051	Ready DNA Bead Stand	

